

## MS4A4A antibody, N-term

**Cat. No. GTX44847**

<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	WB
<b>Reactivity</b>	Human

**Package**  
50 µg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

WB	0.2-2.5 ug/ml
----	---------------

Not tested in other applications.

<b>Calculated MW</b>	25 kDa. ( <a href="#">Note</a> )
----------------------	----------------------------------

## Properties

<b>Form</b>	Liquid
-------------	--------

<b>Buffer</b>	PBS, 2% Sucrose
---------------	-----------------

<b>Preservative</b>	0.09% Sodium azide
---------------------	--------------------

<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
----------------	--

<b>Concentration</b>	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
----------------------	--

<b>Immunogen</b>	A synthetic peptide corresponding to a N-terminal region of Human MS4A4A
------------------	--

<b>Purification</b>	Affinity Purified
---------------------	-------------------

<b>Conjugation</b>	Unconjugated
--------------------	--------------

**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES

**GTX44847 WB Image**

WB analysis of human B lymphoblast cells using GTX44847 MS4A4A antibody at 0.2-1 µg/ml.



For full product information, images and publications, please visit our [website](#).