

C14ORF104 antibody, N-term

Cat. No. GTX44972

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.2-2.5 ug/ml
ICC/IF	Assay dependent

Not tested in other applications.

Calculated MW 91 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a N-terminal region of Human C14ORF104
Purification	Affinity Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX44972 WB Image

WB analysis of HeLa cells using GTX44972 C14ORF104 antibody at 0.2-1 µg/ml.



GTX44972 ICC/IF Image

ICC/IF analysis of human nasal epithelial cells using GTX44972 C14ORF104 antibody at 1:500. Cells were fixed with 4% PFA, permeabilized with 0.2% Triton X-100, blocked with 1% skim milk, and incubated with primary antibody at 4 degree C for overnight.

A : acetylated alpha-tubulin

B : C14ORF104

C : DAPI

D : merge



For full product information, images and publications, please visit our [website](#).