

# DAP12 antibody

**Cat. No. GTX45098**

|                     |            |
|---------------------|------------|
| <b>Host</b>         | Rabbit     |
| <b>Clonality</b>    | Polyclonal |
| <b>Isotype</b>      | IgG        |
| <b>Applications</b> | WB         |
| <b>Reactivity</b>   | Rat        |

**Package**  
100 µg

## Applications

### Application Note

For WB: Use at a dilution of 1:1,000. Optimal dilutions/concentrations should be determined by the end user.

|                      |                                                        |
|----------------------|--------------------------------------------------------|
| <b>Calculated MW</b> | 12 kDa. ( <a href="#">Note</a> )                       |
| <b>Product Note</b>  | DAP12 detected in tubular interstitial nephritis model |

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS                                                                                                                                                                                                                        |
| <b>Preservative</b>  | 0.1% Sodium azide                                                                                                                                                                                                          |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| <b>Immunogen</b>     | Recombinant protein - amino acids 28-40 of rat DAP12.                                                                                                                                                                      |
| <b>Purification</b>  | Protein G affinity purified<br>From polyclonal serum                                                                                                                                                                       |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTx45098 WB Image

Rat DAP12 detected on Day 3, 8, and 10 by Western blot in a tubular interstitial nephritis model



For full product information, images and publications, please visit our [website](#).