

MEK1 antibody

Cat. No. GTX45127

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

For WB: Use at a dilution of 1:2,000. For IP: Use at a dilution of 1:2,000. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW	43 kDa. (Note)
Product Note	MEK1 in RAW cell and RBL transfected with NIK

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein - full length mouse MEK1.
Purification	Protein A affinity purified From polyclonal serum
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX45127 WB Image

Western blot detection of Mek1 protein in RAW cell (1) and RBL transfected with NIK(2).



For full product information, images and publications, please visit our [website](#).