

SUN1 antibody, N-term

Cat. No. GTX45959

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Dog

References (1)

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	2-10 ug/ml

Not tested in other applications.

Calculated MW 90 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a N-terminal region of Human UNC84A
Purification	Affinity Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX45959 WB Image

WB analysis of human fetal lung tissue using GTX45959 UNC84A antibody at 1µg/ml.




GTX45959 IHC-P Image

ICC/IF analysis of C2C12 cells using GTX45959 UNC84A antibody at 1:500.


Green : UNC84A

Blue : DAPI



GTX45959 WB Image

WB analysis of human placenta tissue using GTX45959 UNC84A antibody at 1µg/ml.



GTX45959 WB Image

WB analysis of Jurkat cells using GTX45959 UNC84A antibody at 0.2-1µg/ml.



For full product information, images and publications, please visit our [website](https://www.genetex.com).