

WBP11 antibody, N-term

Cat. No. GTX46466

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IP
Reactivity	Human, Mouse

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.2-2.5 ug/ml
IHC-P	2-10 ug/ml
IP	Assay dependent

Not tested in other applications.

Calculated MW 70 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a N-terminal region of Human WBP11
Purification	Affinity Purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX46466 IHC-P Image


IHC-P analysis of human pineal tissue using GTX46466 WBP11 antibody at 1:100. Left to right : DAPI, WBP11, Merge. Low pH, heat-induced antigen retrieval method utilizing Sodium Citrate buffer was performed.


GTX46466 WB Image

WB analysis of HeLa cells using GTX46466 WBP11 antibody at 0.2-1 µg/ml.


GTX46466 IP Image

IP analysis of mouse brain tissue using GTX46466 WBP11 antibody at 1:500. IP antibody: 6 µg.


GTX46466 WB Image

WB analysis of human NT-2 cells and mouse brain tissue extract using GTX46466 WBP11 antibody at 2 µg/ml.

Lane 1 : human NT-2 cells

Lane 2 : mouse brain tissue



For full product information, images and publications, please visit our [website](https://www.genetex.com).