

# CRELD1 antibody, C-term

**Cat. No. GTX46585**

|                     |            |
|---------------------|------------|
| <b>Host</b>         | Rabbit     |
| <b>Clonality</b>    | Polyclonal |
| <b>Isotype</b>      | IgG        |
| <b>Applications</b> | WB, IHC-P  |
| <b>Reactivity</b>   | Human      |

**Package**  
100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | Assay dependent      |
| IHC-P              | 2-10 ug/ml           |

Not tested in other applications.

**Calculated MW** 45 kDa. ( [Note](#) )

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS, 2% Sucrose                                                                                                                                                                                                            |
| <b>Preservative</b>  | 0.09% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                               |
| <b>Immunogen</b>     | A synthetic peptide corresponding to a C-terminal region of Human CRELD1                                                                                                                                                   |
| <b>Purification</b>  | Protein A purified                                                                                                                                                                                                         |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTX46585 WB Image

WB analysis of HepG2 cells using GTX46585 CRELD1 antibody at 5.0µg/ml.

Lane A : marker

Lane B : HepG2 cells



### GTX46585 IHC-P Image

IHC-P analysis of human kidney tissue using GTX46585 CRELD1 antibody at 4.0-8.0µg/ml.



For full product information, images and publications, please visit our [website](#).