

DPP6 antibody, Internal

Cat. No. GTX46610

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.2-2.5 ug/ml

Not tested in other applications.

Calculated MW 98 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to an Internal region of Human DPP6
Purification	Affinity Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX46610 WB Image

WB analysis of human fetal muscle tissue using GTX46610 DPP6 antibody at 0.2-1µg/ml.



GTX46610 WB Image

WB analysis of mouse brain (from WT mice and DPP6 knockout mice) tissue using GTX46610 DPP6 antibody at 1:1000. (Loading: 20µg)

Lane 1: WT mouse brain extract

Lane 2: DPP6 -/- mouse brain extract

Lane 3: WT mouse brain extract

Lane 4: DPP6 -/- mouse brain extract



For full product information, images and publications, please visit our [website](https://www.genetex.com).