

# MCTP1 antibody, Internal

**Cat. No. GTX46724**

|                     |            |
|---------------------|------------|
| <b>Host</b>         | Rabbit     |
| <b>Clonality</b>    | Polyclonal |
| <b>Isotype</b>      | IgG        |
| <b>Applications</b> | WB, IHC-P  |
| <b>Reactivity</b>   | Human      |

**Package**  
50 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 0.2-2.5 ug/ml        |
| IHC-P              | 2-10 ug/ml           |

Not tested in other applications.

**Calculated MW** 112 kDa. ([Note](#))

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS, 2% Sucrose                                                                                                                                                                                                            |
| <b>Preservative</b>  | 0.09% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                               |
| <b>Immunogen</b>     | A synthetic peptide corresponding to an Internal region of Human MCTP1                                                                                                                                                     |
| <b>Purification</b>  | Affinity Purified                                                                                                                                                                                                          |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTX46724 WB Image

WB analysis of HeLa cells using GTX46724 MCTP1 antibody at 1 µg/ml.



### GTX46724 IHC-P Image

IHC-P analysis of human spleen tissue using GTX46724 MCTP1 antibody at 5 µg/ml.



### GTX46724 IHC-P Image

IHC-P analysis of human spleen tissue using GTX46724 MCTP1 antibody at 5 µg/ml.



For full product information, images and publications, please visit our [website](https://www.genetex.com).