

# SLC36A3 antibody, N-term

**Cat. No. GTX46758**

|                     |            |
|---------------------|------------|
| <b>Host</b>         | Rabbit     |
| <b>Clonality</b>    | Polyclonal |
| <b>Isotype</b>      | IgG        |
| <b>Applications</b> | WB, IHC-P  |
| <b>Reactivity</b>   | Human      |

**Package**  
100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | Assay dependent      |
| IHC-P              | 2-10 ug/ml           |

Not tested in other applications.

**Calculated MW** 52 kDa. ( [Note](#) )

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS, 2% Sucrose                                                                                                                                                                                                            |
| <b>Preservative</b>  | 0.09% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                               |
| <b>Immunogen</b>     | A synthetic peptide corresponding to a N-terminal region of Human SLC36A3                                                                                                                                                  |
| <b>Purification</b>  | Protein A purified                                                                                                                                                                                                         |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |


### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTx46758 IHC-P Image**

IHC-P analysis of human kidney tissue using GTx46758 SLC36A3 antibody at 4.0-8.0µg/ml.


**GTx46758 WB Image**

WB analysis of Jurkat cells using GTx46758 SLC36A3 antibody at 5.0µg/ml.

Lane A : marker

Lane B : Jurkat cells



For full product information, images and publications, please visit our [website](https://www.genetex.com).