

RNF128 antibody, C-term

Cat. No. GTX46906

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.2-2.5 ug/ml

Not tested in other applications.

Calculated MW 47 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a C-terminal region of Human RNF128
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX46906 WB Image

WB analysis of mouse CD4 T cells using GTX46906 RNF128 antibody at 1:2400.

Lane:

- 1: CD4 T cell lysate from WT mice
- 2: CD4 T cell lysate from WT mice treated with anti-CD3/CD28
- 3: CD4 T cell lysate from WT mice treated with FK506
- 4: CD4 T cell lysate from WT mice treated with anti-CD3/CD28 and FK506
- 5: CD4 T cell lysate from TRPV1 KO mice
- 6: CD4 T cell lysate from TRPV1 KO mice treated with anti-CD3/CD28
- 7: CD4 T cell lysate from TRPV1 KO mice treated with FK506
- 8: CD4 T cell lysate from TRPV1 KO mice treated with anti-CD3/CD28 and FK506.



GTX46906 WB Image

WB analysis of HepG2 cells using GTX46906 RNF128 antibody at 1.25µg/ml.



For full product information, images and publications, please visit our [website](https://www.genetex.com).