

SLC19A3 antibody, Internal

Cat. No. GTX46980

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.2-2.5 ug/ml

Not tested in other applications.

Calculated MW 56 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to an Internal region of Human SLC19A3
Purification	Affinity Purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX46980 WB Image


WB analysis of human fetal heart tissue using GTX46980 SLC19A3 antibody at 1µg/ml.


GTX46980 WB Image

WB analysis of DU145 cells using GTX46980 SLC19A3 antibody at 0.2-1µg/ml.


GTX46980 WB Image

WB analysis of human fetal liver tissue using GTX46980 SLC19A3 antibody at 1µg/ml.


GTX46980 WB Image

WB analysis of human ovary tumor tissue using GTX46980 SLC19A3 antibody at 1µg/ml.



For full product information, images and publications, please visit our [website](https://www.genetex.com).