

EIF3G antibody, Internal

Cat. No. GTX47377

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.2-2.5 ug/ml
ICC/IF	Assay dependent

Not tested in other applications.

Calculated MW 36 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to an Internal region of Human EIF3G
Purification	Affinity Purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX47377 WB Image


WB analysis of MCF-7 cells using GTX47377 EIF3G antibody at 1.0µg/ml.


GTX47377 ICC/IF Image


ICC/IF analysis of human brain stem cells using GTX47377 EIF3G antibody at 1:500.

Red: EIF3G

Blue: DAPI


GTX47377 WB Image

WB analysis of 293T cells using GTX47377 EIF3G antibody at 1.0µg/ml.


GTX47377 WB Image

WB analysis of Jurkat cells using GTX47377 EIF3G antibody at 0.2-1µg/ml.



For full product information, images and publications, please visit our [website](https://www.genetex.com).