

KCTD11 antibody, N-term

Cat. No. GTX47451

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC
Reactivity	Human

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.2-2.5 ug/ml
IHC	Assay dependent

Not tested in other applications.

Calculated MW 26 kDa. ([Note](#))

Product Note KO/KD validation is based on published data (PMID: 28465479).

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a N-terminal region of Human KCTD11
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX47451 WB Image

WB analysis of HepG2 cells using GTX47451 KCTD11 antibody at 1.25µg/ml.




GTX47451 WB Image

The data was published in the journal Oncotarget in 2017. [PMID: 28465479](https://pubmed.ncbi.nlm.nih.gov/28465479/)



GTX47451 WB Image

The data was published in the journal Oncotarget in 2017. [PMID: 28465479](https://pubmed.ncbi.nlm.nih.gov/28465479/)



GTX47451 IHC Image

The data was published in the journal Oncotarget in 2017. [PMID: 28465479](https://pubmed.ncbi.nlm.nih.gov/28465479/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).