

## ATF5 blocking peptide

Cat. No. GTX47605-PEP

|         |       |
|---------|-------|
| Species | Human |
|---------|-------|

|         |        |
|---------|--------|
| Package | 100 µg |
|---------|--------|

## Applications

|              |                                                     |
|--------------|-----------------------------------------------------|
| Product Note | The peptide was used in the production of GTX47605. |
|--------------|-----------------------------------------------------|

## Properties

|      |                    |
|------|--------------------|
| Form | Lyophilized powder |
|------|--------------------|

|        |                                                                                                                                         |
|--------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Buffer | Reconstitute with 200µl distilled water to obtain a 0.5mg/ml peptide solution. Lyophilized from 5% (v/v) acetonitrile/H <sub>2</sub> O. |
|--------|-----------------------------------------------------------------------------------------------------------------------------------------|

|              |                  |
|--------------|------------------|
| Preservative | No preservatives |
|--------------|------------------|

|         |                                                                                                                                                                                                                            |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|                 |               |
|-----------------|---------------|
| Region/Sequence | EVYKARSQRTRSC |
|-----------------|---------------|

|      |                                                                                                                                                   |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Note | For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------|

|  |                                                                                                                                                                             |
|--|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product. |
|--|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|



For full product information, images and publications, please visit our [website](#).