

KCNK9 antibody, N-term

Cat. No. GTX47674

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Dog

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.2-2.5 ug/ml
ICC/IF	Assay dependent
IHC-P	2-10 ug/ml

Not tested in other applications.

Calculated MW 42 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a N-terminal region of Human KCNK9
Purification	Affinity Purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX47674 WB Image

WB analysis of HepG2 cells using GTX47674 KCNK9 antibody at 1 µg/ml.



GTX47674 IHC-P Image

IHC-P analysis of human prostate tissue using GTX47674 KCNK9 antibody at 5.0 µg/ml.




GTX47674 ICC/IF Image

ICC/IF analysis of MDCK cells using GTX47674 KCNK9 antibody at 1:100. Cells fixed with 4% PFA for 10 mins, and then blocked with 5% horse serum and 0.1% Triton X-100 for 30 mins.

Green : KCNK9

Blue : DAPI



GTX47674 WB Image

WB analysis of 293T cells using GTX47674 KCNK9 antibody at 1 µg/ml.



For full product information, images and publications, please visit our [website](https://www.genetex.com).