

TLR1 antibody

Cat. No. GTX47794

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IP, ELISA
Reactivity	Mouse, Rat

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500
IHC-P	Assay dependent
IP	1:200
ELISA	1:20000

Not tested in other applications.

Calculated MW 91 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Tris/Glycine, 0.5% BSA, 30% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.55-0.67 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide taken within amino acid region 200-250 on mouse toll-like receptor 1 protein.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx47794 IHC-P Image


IHC-P analysis of rat spleen tissue using GTx47794 TLR1 antibody.

Dilution : 1:100


GTx47794 WB Image

WB analysis of TLR1 recombinant protein using GTx47794 TLR1 antibody.

Dilution : 1:500


GTx47794 IHC-P Image

IHC-P analysis of rat spleen tissue using GTx47794 TLR1 antibody.

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).