

FBXO32 antibody

Cat. No. GTX47819

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IP, ELISA, IHC
Reactivity	Mouse, Rat

References (2)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500
IHC-P	Assay dependent
IP	1:200
ELISA	1:10000
IHC	1:50-1:150

Not tested in other applications.

Calculated MW 42 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Tris/Glycine, 0.5% BSA, 30% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.64-0.72 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide taken within amino acid region 240-280 on Atr1 protein.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTX47819 WB Image

WB analysis of Fbx32 recombinant protein using GTX47819 FBXO32 antibody.


Dilution : 1:500



GTX47819 IHC-P Image

IHC-P analysis of mouse skeletal muscle tissue using GTX47819 FBXO32 antibody.

Dilution : 1:50



GTX47819 IHC-P Image

IHC-P analysis of mouse skeletal muscle tissue using GTX47819 FBXO32 antibody.

Dilution : 1:50



For full product information, images and publications, please visit our [website](https://www.genetex.com).