

Thickvein antibody

Cat. No. GTX47856

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP, ELISA
Reactivity	Drosophila

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500
IP	1:200
ELISA	1:10000

Not tested in other applications.

Properties

Form	Liquid
Buffer	Tris/Glycine, 0.5% BSA, 30% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.7 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide corresponding to unique amino acid sequences on C-terminus of pan Thickvein protein. Immunogenic peptide was covalently modified to achieve the desired antigenicity. Sequence : CVK GFR PPI PSR WQE DD
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note
For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

**GTX47856 WB Image**

WB analysis of Drosophila extracts using GTX47856 Thickvein antibody.

Dilution : 1:500



For full product information, images and publications, please visit our [website](#).