

STRA6 antibody

Cat. No. GTX47858

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IP, ELISA
Reactivity	Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500
IHC-P	Assay dependent
IP	1:200
ELISA	1:10000

Not tested in other applications.

Calculated MW 74 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Tris/Glycine, 0.5% BSA, 30% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.69 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide taken within third large intra cytoplasmic loop within amino acid region 200-300 on rat Stra6 protein.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX47858 IHC-P Image

IHC-P analysis of mouse kidney tissue using GTX47858 STRA6 antibody.
Dilution : 1:100




GTX47858 IHC-P Image

IHC-P analysis of mouse kidney tissue using GTX47858 STRA6 antibody.
Dilution : 1:100



GTX47858 IHC-P Image

IHC-P analysis of mouse spleen tissue using GTX47858 STRA6 antibody.
Dilution : 1:100



GTX47858 IHC-P Image

IHC-P analysis of mouse kidney tissue using GTX47858 STRA6 antibody.
Dilution : 1:100



GTX47858 IHC-P Image

IHC-P analysis of mouse spleen tissue using GTX47858 STRA6 antibody.
Dilution : 1:100



For full product information, images and publications, please visit our [website](#).