

SLC7A10 blocking peptide

Cat. No. GTX47874-PEP

Application	Neutralizing/Inhibition
Species	Rat

Package
250 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Neutralizing/Inhibition	Assay dependent

Not tested in other applications.

PROPERTIES

Form	Liquid
Buffer	PBS, 10% DMSO
Preservative	No preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	2.5 mg/ml (Please refer to the vial label for the specific concentration.)
Region/Sequence	Synthetic peptide taken within amino acid region 1-50 on rat Asc-type amino acid transporter 1 protein.
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).