

Aquaporin 12 antibody

Cat. No. GTX47925

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IP, ELISA
Reactivity	Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500
ICC/IF	1:50-1:100
IHC-P	Assay dependent
IP	1:200
ELISA	1:10000

Not tested in other applications.

Calculated MW 31 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Tris/Glycine, 0.5% BSA, 30% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.61-0.67 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide taken within amino acid region 240-290 on mouse Aquaporin 12 protein.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTX47925 WB Image

WB analysis of Aquaporin 12 recombinant protein using GTX47925 Aquaporin 12 antibody.


Dilution : 1:500



GTX47925 IHC-P Image

IHC-P analysis of mouse kidney tissue using GTX47925 Aquaporin 12 antibody.

Dilution : 1:100



GTX47925 IHC-P Image

IHC-P analysis of rat kidney tissue using GTX47925 Aquaporin 12 antibody.

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).