

Caspase Family Inhibitor, Z-VAD(OMe)-FMK

Cat. No. GTX47941
Package
20 µl

PRODUCT
Summary Cell permeable peptide inhibitor for Caspase family

Properties
Form Liquid

Buffer 10mM DMSO

Preservative No preservative

Storage Store at -20°C. Protect from light and moisture.

Concentration 10 mM (Please refer to the vial label for the specific concentration.)

Purity Single spot by TLC

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2025 / 06 / 13 Page 1 of 1