

Caspase-9 Inhibitor, Z-LE(OMe)HD(OMe)-FMK

Cat. No. GTX47945

Package
20 µl

PRODUCT

Summary Cell permeable peptide inhibitor for Caspase 9

Properties

Form Liquid

Buffer 10mM DMSO

Preservative No preservatives

Storage Store at -20°C. Protect from light and moisture.

Concentration 10 mM (Please refer to the vial label for the specific concentration.)

Purity >98% by HPLC

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).