

Human IL12 protein

Cat. No. GTX47956-pro

Species	Human
----------------	-------

Package	10 µg
----------------	-------

Properties

Form	Lyophilized powder
-------------	--------------------

Buffer	Batch dependent (please contact us for details)
---------------	---

Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
----------------	---

Region/Sequence	p40 subunit: IWELKKDVYV VELDWPDPAP GEMVLTCDT PEEDGITWTL DQSSEVLGSG KTLTIQVKEF GDAGQYTCHK GGEVLSHSLLL LHKKEDGIW STDILKDQKE PKNKTLRCE AKNYSGRFTC WWLTTISTDL TFSVKSSRGS SDPQGVTCGA ATLSAERVRG DNKEYEYSVE CQEDSACPAA EESLPIEMV DAVHKLKYEN YTSSFFIRDI IKPDPPKNLQ LKPLKNSRQV EVSWEYPDTW STPHSYFSLT FCVQVQGKSK REKKDRVFTD KTSATVICRK NASISVRAQD RYSSSWSEW ASVPCS
------------------------	---

Expression System	Baculovirus (High-5 cells)
--------------------------	----------------------------

Purity	> 98% by SDS-PAGE and HPLC.
---------------	-----------------------------

Endotoxin	< 1 EU/µg
------------------	-----------

Conjugation	Unconjugated
--------------------	--------------

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
-------------	---

Note	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
-------------	---



For full product information, images and publications, please visit our [website](#).