

Rat Soluble RANK Ligand protein (active)

Cat. No. GTX47970-pro

Application Functional Assay**Species** Rat**Package**

10 µg

APPLICATION

Application Note

Biological activity was determined by its ability to induce osteoclast formation in RAW264.7 cells using a concentration range of 5.0-10.0 ng/ml.

PROPERTIES

Form Lyophilized powder**Buffer** Batch dependent (please contact us for details)**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** PAMMEGSWLD VARRGKPEAQ PFAHLTINAA DIPSGSHKVS LSSWYHDRGW AKISNMTLSN GKLRVNQDGF YYLYANICFR
HHETSGSVPA DYQLMVVWV KTSIKIPSSH NLMKGGSTKN WSGNSEFHFY SINVGFFKL RAGEEISVQV SNPSLLDPDQ
DATYFGAFKV QDID**Expression System** E. coli**Purity** > 98% by SDS-PAGE and HPLC.**Endotoxin** < 1 EU/µg**Conjugation** Unconjugated**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).