

Mouse Exodus-2 protein (active)

Cat. No. GTX47986-pro**Applications** Functional Assay**Species** Mouse**Package**

20 µg

Applications

Application Note

Biological activity was determined by its ability to chemoattract total murine T cell population using a concentration range of 10.0-100.0 ng/ml.

Properties

Form Lyophilized powder**Buffer** Batch dependent (please contact us for details)**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** SDGGGQDCCL KYSQKKIPYS IVRGYRKQEP SLGCPIAIL FLPRKHSKPE LKANPEEGWV QNLMRRLDQP PAPGKQSPGC
RKNRGTSKSG KKGKGSKGCR KTEQTQPSRG**Expression System** E. coli**Purity** > 98% by SDS-PAGE and HPLC.**Endotoxin** < 1 EU/µg**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).