

Human Neuritin protein (active)

Cat. No. GTX48012-pro

Applications Functional Assay

Species Human

Package

20 µg

Applications

Application Note

Biological activity was determined by a cell proliferation using rat C6 cells. The ED₅₀ for this effect is 20.0-30.0 ng/ml.

Properties

Form Lyophilized powder

Buffer Batch dependent (please contact us for details)

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence MAGKCDAVFK GFSDCLLKG DSMANYPQGL DDKTNIKTVC TYWEDFHSCT VTALDTCQEG AKDMWDKLRK ESKNLNIQGS LFELCGSGN

Expression System E. coli

Purity > 98% by SDS-PAGE and HPLC.

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).