

Human Neuritin protein (active)

Cat. No. GTX48012-pro**Applications** Functional Assay**Species** Human**Package**

20 µg

Applications

Application Note

Biological activity was determined by a cell proliferation using rat C6 cells. The ED₅₀ for this effect is 20.0-30.0 ng/ml.

Properties

Form Lyophilized powder**Buffer** Batch dependent (please contact us for details)**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** MAGKCDAVFK GFSDCLLKG DSMANYPQGL DDKTNIKTVC TYWEDFHSCT VTALTDCEG AKDMWDKLRK ESKNLNIQGS
LFELCGSGN**Expression System** E. coli**Purity** > 98% by SDS-PAGE and HPLC.**Endotoxin** < 1 EU/µg**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).