

Human TACI protein (active)

Cat. No. GTX48045-pro

Applications Functional Assay**Species** Human**Package**

20 µg

Applications

Application Note

Biological activity was determined by its ability to block human BAFF induced T2B cell survival using a concentration range of 1.0-3.0 µg/ml.

Properties

Form Lyophilized powder**Buffer** Batch dependent (please contact us for details)**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** MSGLGSRRRG GRSRVDQEER FPQGLWTGVA MRSCPPEQYW DPLLGTCTMSC KTICNHQSQR TCAAFCSRSL CRKEQGKFYD
HLLRDCISCA SICGQHPKQC AYFCENKLRS PVNLPELRR QRSGEVENNS DNSGRYQGLE HRGSEASPAL PGLKLSADQV**Expression System** E. coli**Purity** > 98% by SDS-PAGE and HPLC.**Endotoxin** < 1 EU/µg**Conjugation** Unconjugated**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).