

# Human DKK1 protein

**Cat. No. GTX48048-pro**

**Species** Human

**Package**  
10 µg

## Properties

**Form** Lyophilized powder

**Buffer** Batch dependent (please contact us for details)

**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

**Region/Sequence** TLNSVLNSNA IKNLPPPLGG AAGHPGSAVS AAPGILYPGG NKYQTIDNYQ PYPCAEDEEC GTDEYCASPT RGGDAGVQIC  
LACRKRKRRC MRHAMCCPGN YCKNGICVSS DQNHFRGEIE ETITESFGND HSTLDGYSRR TTLSSKMYHT KGQEGSVCLR  
SSDCASGLCC ARHFWKICK PVLKEGQVCT KHRKKGSHGL EIFQRCYCGE GLSCRIQKDH HQASNSSRLH TCQRH

**Expression System** HEK293 cells

**Purity** > 97% by SDS-PAGE and HPLC.

**Endotoxin** < 1 EU/µg

**Conjugation** Unconjugated

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).