

Human DKK1 protein

Cat. No. GTX48048-pro

Species Human

Package
10 µg

Properties

Form Lyophilized powder

Buffer Batch dependent (please contact us for details)

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence TLNSVLNSNA IKNLPPPLGG AAGHPGSAVS AAPGILYPGG NKYQTIDNYQ PYPCAEDEEC GTDEYCASPT RGGDAGVQIC
LACRRRRKRC MRHAMCCPGN YCKNGICVSS DQNHFRGEIE ETITESFGND HSTLDGYSRR TTLSSKMYHT KGQEGSVCLR
SSDCASGLCC ARHFWKICK PVLKEGQVCT KHRKGGSHGL EIFQRCYCGE GLSCRIQKDH HQASNSSRLH TCQRH

Expression System HEK293 cells

Purity > 97% by SDS-PAGE and HPLC.

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).