

# Mouse CTACK protein (active)

**Cat. No. GTX48055-pro**

<b>Applications</b>	Functional Assay
<b>Species</b>	Mouse

**Package**  
20 µg

## Applications

### Application Note

Biological activity was determined by its ability to chemoattract CXCR3 transfected HEK/293 cells using a concentration range of 10.0-100.0 ng/ml.

## Properties

<b>Form</b>	Lyophilized powder
<b>Buffer</b>	Batch dependent (please contact us for details)
<b>Storage</b>	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
<b>Region/Sequence</b>	LPLPSSTSCC TQLYRQPLPS RLLRRIVHME LQEADGDCHL QAWVLHLARR SVCVHPQNRS LARWLERQ GK RLQGTVP SLN LVLQKKMYSN PQQQN
<b>Expression System</b>	E. coli
<b>Purity</b>	> 98% by SDS-PAGE and HPLC.
<b>Endotoxin</b>	< 1 EU/µg
<b>Conjugation</b>	Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).