

Mouse IL11 protein (active)

Cat. No. GTX48066-pro

Applications Functional Assay

Species Mouse

Package

10 µg

Applications

Application Note

The ED₅₀ was determined by the dose-dependent stimulation of the proliferation of murine T11 cells is ≤ 2.0 ng/ml, corresponding to a specific activity of $\geq 5 \times 10^5$ units/mg.

Properties

Form Lyophilized powder

Buffer Batch dependent (please contact us for details)

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence MPGPPAGSPR VSSDPRADLD SAVLLTRSL ADTRQLAAQM RDKFPADGDH SLDSLPTLAM SAGTLGSLQL PGVLRRLRVD
LMSYLRHVQW LRRAGGPSLK TLEPELGALQ ARLERLLRRL QLLMSRLALP QAAPDQPVIP LGPPASAWGS IRAAHAILGG
LHLTDWAVR GLLLLKTRL

Expression System E. coli

Purity > 98% by SDS-PAGE and HPLC.

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).