

Mouse Interferon Gamma protein (active)

Cat. No. GTX48067-pro

Applications Functional Assay**Species** Mouse**Package**

100 µg

Applications

Application Note

The ED₅₀ was determined by a cytopathic affect inhibition assay with murine NIH 3T3 cells challenged with VSV virus was 0.1 ng/ml, corresponding to a specific activity of 1 x 10⁷ units/mg.

Properties

Form Lyophilized powder**Buffer** Batch dependent (please contact us for details)**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** MHGTVIESLE SLNNYFNSSG IDVEEKSFL DIWRNWQKDG DMKILQSQII SFYLRLEVL KDNQAISNNI SVIESHLITT FFSNSKAKKD AFMSIAKFEV NNPQVQRQAF NELIRVVHQL LPESLRKRK RSRC**Expression System** E. coli**Purity** > 98% by SDS-PAGE and HPLC.**Endotoxin** < 1 EU/µg**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).