

Mouse RCAN2 Isoform 2 protein, His tag

Cat. No. GTX48083-pro

Species	Mouse
----------------	-------

Package 50 µg

Properties

Form	Liquid
-------------	--------

Buffer	10mM Tris, 100mM Sodium Phosphate, 500mM NaCl, 25mM Imidazole, 2mM MgCl ₂ , 10% Glycerol
---------------	---

Preservative	No preservatives
---------------------	------------------

Storage	Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
----------------	--

Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
----------------------	--

Region/Sequence	C-terminal His-Tag; a.a. 2-197
------------------------	--------------------------------

Expression System	E. coli
--------------------------	---------

Purity	~95%
---------------	------

Conjugation	Unconjugated
--------------------	--------------

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
-------------	---



For full product information, images and publications, please visit our [website](#).