

Human DNAL4 protein

Cat. No. GTX48084-pro

Species	Human
----------------	-------

Package 50 µg

PROPERTIES

Form	Lyophilized powder
-------------	--------------------

Buffer	Batch dependent (please contact us for details)
---------------	---

Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
----------------	---

Region/Sequence	a.a.1-105
------------------------	-----------

Expression System	E. coli
--------------------------	---------

Purity	> 90%
---------------	-------

Conjugation	Unconjugated
--------------------	--------------

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
-------------	---

	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
--	---



For full product information, images and publications, please visit our [website](#).