

Human XBP1 protein, His tag

Cat. No. GTX48114-pro

Species Human

Package
50 µg

Properties

Form Liquid

Buffer PBS, 0.1% SDS

Preservative No preservatives

Storage Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration Batch dependent (Please refer to the vial label for the specific concentration.)

Region/Sequence N-terminal His-Tag; a.a. 2-160

Expression System E. coli

Purity ~95%

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).