

Human SGTA protein

Cat. No. GTX48139-pro

| | |
|----------------|-------|
| Species | Human |
|----------------|-------|

| |
|-------------------------|
| Package 50 µg |
|-------------------------|

PROPERTIES

| | |
|-------------|--------------------|
| Form | Lyophilized powder |
|-------------|--------------------|

| | |
|---------------|---|
| Buffer | Batch dependent (please contact us for details) |
|---------------|---|

| | |
|----------------|---|
| Storage | Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles. |
|----------------|---|

| | |
|------------------------|-----------|
| Region/Sequence | a.a.1-313 |
|------------------------|-----------|

| | |
|--------------------------|---------|
| Expression System | E. coli |
|--------------------------|---------|

| | |
|---------------|-----|
| Purity | .95 |
|---------------|-----|

| | |
|--------------------|--------------|
| Conjugation | Unconjugated |
|--------------------|--------------|

| | |
|-------------|---|
| Note | For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. |
|-------------|---|

| | |
|--|---|
| | Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product. |
|--|---|



For full product information, images and publications, please visit our [website](#).