

Human PEX19 protein

Cat. No. GTX48152-pro

Species	Human
----------------	-------

Package
50 µg

Properties

Form	Lyophilized powder
-------------	--------------------

Buffer	Batch dependent (please contact us for details)
---------------	---

Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
----------------	---

Region/Sequence	a.a.1-299
------------------------	-----------

Expression System	E. coli
--------------------------	---------

Purity	.95
---------------	-----

Conjugation	Unconjugated
--------------------	--------------

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
-------------	---

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).