

# Human Neurturin protein (active)

**Cat. No. GTX48169-pro**

<b>Applications</b>	Functional Assay
<b>Species</b>	Human

**Package**  
20 µg

## Applications

### Application Note

Human Neurturin at a concentration of 100 ng/ml will support the survival of 65% of newborn rat sympathetic neurons.

## Properties

<b>Form</b>	Lyophilized powder
<b>Buffer</b>	Batch dependent (please contact us for details)
<b>Storage</b>	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
<b>Region/Sequence</b>	MARLGARPCG LRELEVRVSE LGLGYASDET VLFYRCAGAC EAAARVYDLG LRRLRQRRRL RRERVRAQPC CRPTAYEDEV SFLDAHSRYH TVHELSAREC ACV
<b>Expression System</b>	E. coli
<b>Purity</b>	> 98% by SDS-PAGE and HPLC.
<b>Endotoxin</b>	< 1 EU/µg
<b>Conjugation</b>	Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).