

Human Neurturin protein (active)

Cat. No. GTX48169-pro

Applications Functional Assay

Species Human

Package

20 µg

Applications

Application Note

Human Neurturin at a concentration of 100 ng/ml will support the survival of 65% of newborn rat sympathetic neurons.

Properties

Form Lyophilized powder

Buffer Batch dependent (please contact us for details)

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence MARLGARPCG LRELEVRVSE LGLGYASDET VLFYRCAGAC EAAARVYDLG LRRLRQRRRL RRERVRAQPC CRPTAYEDEV
SFLDAHSRYH TVHELSAREC ACV

Expression System E. coli

Purity > 98% by SDS-PAGE and HPLC.

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).