

Mouse NELL2 protein, His tag

Cat. No. GTX48171-pro

Species	Mouse
----------------	-------

Package	50 µg
----------------	-------

Properties

Form	Liquid
-------------	--------

Buffer	10mM Tris, 100mM Sodium Phosphate, 500mM NaCl, 25mM Imidazole, 2mM MgCl ₂ , 10% Glycerol
---------------	---

Preservative	No preservatives
---------------------	------------------

Storage	Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
----------------	--

Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
----------------------	--

Region/Sequence	C-terminal His-Tag; a.a. 22-238
------------------------	---------------------------------

Expression System	E. coli
--------------------------	---------

Purity	~95%
---------------	------

Conjugation	Unconjugated
--------------------	--------------

Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>
-------------	---



For full product information, images and publications, please visit our [website](#).