

# Human IGFBP7 protein

**Cat. No. GTX48190-pro**

**Species** Human

**Package**

25 µg

## Properties

**Form** Lyophilized powder

**Buffer** Batch dependent (please contact us for details)

**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

**Region/Sequence** MSSSDTCGPC EPASCPPLPP LGCLGETRD ACGCCPMCAR GEGEPCGGGG AGRGYCAPGM ECVKSRKRRK GKAGAAAGGP  
GVSGVCVCKS RYPVCGSDGT TYPGCGQLRA ASQRAESRGE KAITQVSKGT CEQGSPVTP PKDIWNVGTGA QVYLSCEVIG  
IPTPVLWINK VKRGHYGVQR TELLPGDRDN LAIQTRGGPE KHEVTGWVLV SPLSKEDAGE YECHASNSQG QASASAKITV  
VDALHEIPVK KGEAEL

**Expression System** E. coli

**Purity** > 98% by SDS-PAGE and HPLC.

**Endotoxin** < 1 EU/µg

**Conjugation** Unconjugated

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).