

Influenza A virus HA (Hemagglutinin) protein, His and StrepII tag

Cat. No. GTX48233-pro

Species	Influenza A virus
----------------	-------------------

Package	50 µg
----------------	-------

Properties

Form	Liquid
-------------	--------

Buffer	PBS, 0.1% SDS
---------------	---------------

Preservative	0.02% Sodium azide
---------------------	--------------------

Storage	Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
----------------	--

Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
----------------------	--

Region/Sequence	N-terminal His-Tag and StrepII-Tag; a.a. 17-338
------------------------	---

Expression System	E. coli
--------------------------	---------

Purity	~95%
---------------	------

Conjugation	Unconjugated
--------------------	--------------

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).