

# Influenza A virus HA (Hemagglutinin) protein, His and StreptII tag

**Cat. No. GTX48233-pro**

**Species** Influenza A virus

**Package**  
50 µg

## Properties

**Form** Liquid

**Buffer** PBS, 0.1% SDS

**Preservative** 0.02% Sodium azide

**Storage** Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

**Concentration** Batch dependent (Please refer to the vial label for the specific concentration.)

**Region/Sequence** N-terminal His-Tag and StreptII-Tag; a.a. 17-338

**Expression System** E. coli

**Purity** ~95%

**Conjugation** Unconjugated

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).