

Human IL28A protein (active)

Cat. No. GTX48237-pro

Applications	Functional Assay
Species	Human

Package
20 µg

Applications

Application Note

Biological activity was determined by its ability to activate STAT following receptor-ligand interaction.

Properties

Form	Lyophilized powder
Buffer	Batch dependent (please contact us for details)
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
Region/Sequence	PVARLHGALP DARGCHIAQF KSLSPQELQA FKRAKDALEE SLLKDCRCH SRLFPRTWDL RQLQVRERPM ALEAELALT KVLEATADTD PALVDVLDQP LHTLHHILSQ FRACIQPQPT AGPRTRGRLH HWLYRLQEAP KKESPGCLEA SVTFNLFRL TRDLNCSVASG DLCV
Expression System	E. coli
Purity	> 98% by SDS-PAGE and HPLC.
Endotoxin	< 1 EU/µg
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).