

Mouse Soluble RANK Ligand protein (active)

Cat. No. GTX48283-pro

Applications Functional Assay

Species Mouse

Package

10 µg

Applications

Application Note

Biological activity was measured by its ability to induce osteoclast formation on murine RAW264.7 cells using a concentration of 5.0-10.0 ng/ml.

Observed MW (kDa) 19.4 kDa.

Properties

Form Lyophilized powder

Buffer Batch dependent (please contact us for details)

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence PAMMEGSWLD VAQRGKPEAQ PFAHLTINAA SIPSGSHKVT LSSWYHDRGW AKISNMTLSN GKLRVNQDGF YYLYANICFR
HHETSGSVPT DYQLMVYVW KTSIKIPSSH NLMKGGSTKN WSGNSEHFY SINVGGFKL RAGEEISIQV SNPSLLDPDQ
DATYFGAFKV QDID

Expression System E. coli

Purity > 98% by SDS-PAGE and HPLC.

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).