

Rat IL3 Beta protein (active)

Cat. No. GTX48303-pro

Applications Functional Assay**Species** Rat**Package**

20 µg

Applications

Application Note

The ED₅₀ was determined by the dose-dependent stimulation of the proliferation of thymidine uptake by murine MC-9 cells is < 10 ng/ml, corresponding to a specific activity of > 1.0 x 10⁵ units/mg.

Properties

Form Lyophilized powder**Buffer** Batch dependent (please contact us for details)**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** MISDRGSDAH HLLRTLDCRT IALEILVKLP YPQVSGLNNS DDKANLRNST LRRVNLDEFL KSQEEFDSQD TTDIKSKLQK
LKCCIPAAAS DSVLPGVYNK DLDDFKKKLR FYVIHLKDLQ PVSISRPPQP TSSSDNFRPM TVEC**Expression System** E. coli**Purity** > 98% by SDS-PAGE and HPLC.**Endotoxin** < 1 EU/µg**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).