

## Mouse IL17A protein (active)

**Cat. No. GTX48307-pro**

<b>Applications</b>	Functional Assay
<b>Species</b>	Mouse

**Package**  
25 µg

## Applications

**Application Note**

Measured by its ability to induce IL-6 production by NIH 3T3 cells.

## Properties

<b>Form</b>	Lyophilized powder
<b>Buffer</b>	Batch dependent (please contact us for details)
<b>Storage</b>	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
<b>Region/Sequence</b>	MAAIIQSSA CPNTEAKDFL QNVKVNLIKVF NSLGAKVSSR RPSDYLNIRST SPWTLHRNED PDRYPSVIWE AQCRHQRCVN AEGKLDHHMN SVLIQQEILV LKREPESCPF TFRVEKMLVG VGCTCVASIV RQAA
<b>Expression System</b>	E. coli
<b>Purity</b>	> 98% by SDS-PAGE and HPLC.
<b>Endotoxin</b>	< 1 EU/µg
<b>Conjugation</b>	Unconjugated

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).